



FL-10000U

Humidification
Designed for Flexibility

Adjustable Temperature

Tailor temperature settings to fit patient needs with precision.

Simplified Workflow


User-friendly touchscreen with intuitive controls. Reduced reprocessing, when used with the Flexicare single patient use temperature probe.

Standard Connection Design

Minimize workflow disruption with a familiar setup.



Single-patient use disposable temperature probe



Product Specifications

Dimensions	W 135mm x D 170mm x H 156mm
Weight	3.3kg
Display	5.6" Touch Screen TFT LCD
Accuracy	±2°C between 30 to 42°C (in humid air)
Heater Plate Thermal Cutout	115±3°C
Ambient Air Temperature Operating Range	18-26°C
Humidity Output (BTPS)	Invasive: 33mg/l - Non-invasive: 16mg/l
Flow Rate	Invasive: Min 1.5lpm at PEEP 3 cmH ₂ O - Non-Invasive: Max 60lpm

FL-10000U Respiratory Humidifier

SKU	Description	QTY
FL-10000U	FL-10000U Respiratory Humidifier USA 110-120V with NEMA 5-15P 3 pin plug	1

FL-10000U Approved Breathing Circuits

SKU	Description	QTY
038-31-101U	Adult Single Heated Wire Breathing System 22mm 1.6m Kit 1	10
038-31-203U	Adult Dual Heated Wire Breathing System 22mm 1.6m Kit 3	10
038-31-252U	Adult Ventimyst Dual Heated Wire Breathing System 1.6m (63") Kit 1	10
038-31-402U	Adult CPAP /Hot Single Limb Heated Wire Breathing System 22mm 1.6m Kit 2 1.6mm	10
038-31-501U	Adult Shielded Single Limb Heated Wire Breathing System 1.6m (63") Kit 1	10
038-33-203U	Neonatal Dual Heated Wire Breathing System 1.6m Kit 3	10
038-33-252U	Neonatal Ventimyst Dual Heated Wire Breathing System 1.6m (63") Kit 1	10
038-33-402U	Neonatal Single Limb Heated Wire Breathing System 1.6m Kit 1	10
038-31-742U	Double Valve Auto-Fill Humidification Chamber	20

FL-10000U Approved Accessories

SKU	Description	QTY
FL-10000-01U	Single Heated Wire Humidifier Adapter Lead	1
FL-10000-02U	Dual Heated Wire Humidifier Adapter Lead	1
FL-10000-03U	Temperature Probe Lead 1.5m	1
FL9000-04U	Single Use Temperature Probe Lead 1.5m	20
FL2000-01U	Mounting Accessory - Pole Clamp (30-50mm)	1
FL2000-02U	Mounting Accessory - C-Clamp (19-40mm)	1
FL2000-03U	Mounting Accessory - Wall Bracket	1
FL2000-04U	Mounting Accessory - Rail Bracket	1
FL2000-05U	Mounting Accessory - Offset Rail Bracket	1



Flexicare Inc.

25671 Commercentre Drive • Lake Forest • CA 92630 USA

949-450 9999 • inquiries@flexicare.com

www.flexicare.com

MM1000-085A | FC9853 | 052026